

**3TEKS TEKSTİL GIDA İTHALAT İHRACAT SANAYİ VE TİCARET LİMİTED ŞİRKETİ**

Gaziantep (Factory) Tel: 00.90.342 337 24 16 İstanbul (Office) Tel: 00.90.216 384 24 48 Ankara (Office) Tel: 00.90.312 442 15 79  
00.90.342 337 24 19 00.90.216 384 24 49 00.90.312 442 15 89  
Fax: 00.90.342 337 21 11 Fax: 00.90.216 384 24 57

2. Organize Sanayi Bölgesi Muammer Güler Bulvarı No:1 GAZIANTEP /TÜRKİYE  
Web: www.3teks.com.tr Email : info@3teks.com.tr




**DATA SHEET**

**DISPOSABLE STERILE STANDARD SURGICAL GOWNS**

CONTENT	SIZE
Standard Surgical Gowns	S – M – L – XL - XXL



- **Fabric:** SSMMS
- **Fabric Weight:** 45 gr/m<sup>2</sup>
- **Sterilization Type:** EO Sterile
- **Packing:** Flat
- **Fabric for medical intention**
- Bacteria-proof
- Anti-microbial,
- Anti-allergenic
- Anti-static
- Air permeable
- Repellant
- Non-linting
- Easy to apply
- **50 Pcs/Box**
- **Box Dimensions: 40\*60\*44 cm**
- **Lengths for Each Sizes;**
  - Small – 120 cm
  - Medium – 130 cm
  - Large – 140 cm
  - XLarge – 150 cm
  - XXLarge – 160 cm



TSE EN 13795-1  
EN ISO 9001:2000  
0197 EN ISO 13485:2003



**3TEKS TEKSTİL GIDA İTHALAT İHRACAT SANAYİ VE TİCARET LİMİTED ŞİRKETİ**

Gaziantep (Factory) Tel: 00.90.342 337 24 16 İstanbul (Office) Tel: 00.90.216 384 24 48 Ankara (Office) Tel: 00.90.312 442 15 79  
00.90.342 337 24 19 00.90.216 384 24 49 00.90.312 442 15 89  
Fax: 00.90.342 337 21 11 Fax: 00.90.216 384 24 57

2. Organize Sanayi Bölgesi Muammer Güler Bulvarı No:1 GAZIANTEP /TÜRKİYE  
Web: www.3teks.com.tr Email : info@3teks.com.tr




**DATA SHEET**

**DISPOSABLE STERILE REINFORCED SURGICAL GOWNS**

CONTENT	SIZE
Reinforced Surgical Gowns	S – M – L – XL – XXL



- **Fabric:** SSMMS
- **Fabric Weight:** 35 gr/m<sup>2</sup>
- **Reinforcement** (Thermobond PE+NW) on the front side and sleeves
- **Sterilization Type:** EO Sterile
- **Packing:** Flat
- **Fabric for medical intention**
- Bacteria-proof
- Anti-microbial,
- Anti-allergenic
- Anti-static
- Air permeable
- Repellant
- Non-linting
- Easy to apply
- **50 Pcs/Box**
- **Box Dimensions: 40\*60\*44 cm**
- **Lengths for Each Sizes;**
  - Large – 140 cm
  - XLarge – 150 cm
  - XXLarge – 160 cm



TSE EN 13795-1  
EN ISO 9001:2000  
0197 EN ISO 13485:2003

