




PRODUCT										QUANTITY			MEASURES	PALLET			
Product Images	Brand	Description	EC Mark	Available Measurers				Code	EAN Code	Nr. pcs per bag	Nr. bag per carton	Nr. pcs per carton	Carton in mm	Layers	Cartons		
<b>INFUSION THERAPY DEVICES</b>																	
<i>Sets for infusional therapies in PVC</i>																	
	<b>INFUSION-SET</b>	Standard Infusion set in PVC, composed by standard or Z-Type ABS spike with air vent, in PVC, with 0.45 micron filter and polyethylene cap; 6 cm of chamber in PVC without above cited phthalates or PVC, and with micro drip (60 drops=1±0,1 ml) or macro drip (20 drops=1±0,1 ml) and with 15 micron filter solution; the tube, in PVC, is memory-less, with external and inner diameter about 4,1x3 mm; ABS roller; isoprene and latex free injection point; ABS Connector with PE cap. Possibility of inserting ABS flow regulator drum shape and of mounting 21G Needle (0,8x38 mm) with its PE protection cap, on the luer lock connector.	In compliance with Directive 93/42EEC according to Medical Device as modified by 2007/47/EEC	Tube	Air vent	Injection site	Roller	Connector	Blister (B) Polibag (P)								
				150cm	Standard	Y	Standard	Luer	P	DEF2YL	8032764150372	25	20	500	580x380x420	4	16
				150cm	Standard	Y	Standard	Luer	B	DEF2YLB	8032764159399	20	18	360	580x380x380	5	20
				150cm	Standard	Y	Standard	Luer Lock	P	DEF2YLL	8032764150358	25	20	500	580x380x420	4	16
				150cm	Standard	Y	Standard	Luer Lock	B	DEF2YLLB	8032764151201	20	18	360	580x380x380	5	20
				150cm	Standard	Y	Standard	Luer Lock-21G	B	DEF2YLL21B	8034055918173	20	18	360	580x380x380	5	20
				200cm	Standard	Y	Standard	Luer Lock	P	DEF2YLL2	8032764154974	25	20	500	580x380x420	4	16
				200cm	Standard	Y	Standard	Luer Lock	B	DEF2YLL2B	8034055913192	20	18	360	580x380x380	5	20
				150cm	Standard	Barilotto	Standard	Luer	P	DEF2BL	8032764150389	25	20	500	580x380x420	4	16
				150cm	Standard	Barilotto	Standard	Luer 21G	P	DEF2BL21	8032764151133	25	20	500	580x380x420	4	16
				150cm	Standard	Barilotto	Standard	Luer	B	DEF2BLB	8034055914373	20	18	360	580x380x380	5	20
				150cm	Standard	Barilotto	Standard	Luer Lock	P	DEF2BLL	8032764150365	25	20	500	580x380x420	4	16
				150cm	Standard	Barilotto	Standard	Luer Lock	B	DEF2BLLB	8032764151973	20	18	360	580x380x380	5	20
				150cm	Standard	Barilotto	Standard	Luer Lock-21G	P	DEF2BLL21	8032764150631	25	20	500	580x380x420	4	16
				150cm	Standard	Barilotto	Standard	Luer Lock-21G	B	DEF2BLL21B	-	20	18	360	580x380x380	5	20
				150cm	Standard	Barilotto	High Precision	Luer Lock	B	DEF2BRLLB	8034055922347	20	18	360	580x380x380	5	20
				150cm	Standard	Barilotto	Ad alta precisione con micro gocciolatore	Luer Lock	B	DEM2BRLLB	8034055922330	20	18	360	580x380x380	5	20
150cm LIGHT SENSITIVE	Standard	Barilotto	Standard	Luer Lock	P	DEF2BLLF	-	25	20	500	580x380x420	4	16				
200cm	Standard	Barilotto	Standard	Luer Lock	P	DEF2BLL2	8032764154967	25	20	500	580x380x420	4	16				
200cm	Standard	Barilotto	Standard	Luer Lock	B	DEF2BLL2B	8032764159405	20	18	360	580x380x380	5	20				