

Betalancet Single-Use-Safety Lancets⁹⁾

- Extraction of capillary blood by contact-activation
- Available size: 28G, 1.8 mm (other sizes on request)
- Automatic Interlock immediately after pricking
- User Protection against stab wounds
- Easy Use
- Available packaging unit: 100 pcs.



* Listed by BGW (Professional Association for Health Service and Welfare Work) as a safe product.

** TRBA 250: Technical Rules for Biological Agents - Biological Agents in Healthcare and Welfare Work, GMBI 2014, Nr. 10/11, March 2014



⁹⁾ Trademark owned and distributed by IME-DC GmbH