

OctaCare

Wound Care



OctaCare Transparent PU Surgical Tape

Used for fixation of medical devices on the skin. The carrier is thin and transparent. Hypoallergenic.



- ◆ Waterproof due to its polyurethane structure
- ◆ Provides comfort of having shower
- ◆ Air permeable
- ◆ Preventing bacterial contamination
- ◆ Can be easily sized by the scissors
- ◆ The flexible structure does not restrict the daily movement of the patient
- ◆ Easy fixation
- ◆ X ray permeable

Code	Size	Quantity / Box	Quantity / Carton
14510	5m x 10cm	1 pcs	60 pcs
14110	10m x 10cm	1 pcs	60 pcs
14115	10m x 15cm	1 pcs	40 pcs

Store in a cool and dry place.

OctaCare Transparent Wound Dressing

Used on all types of wounds. Ideal for protecting post operative wound area and sutured skin. Preserving the wound area from water. Provides comfort of having shower. Ideal for people with sensitive skin. Suitable for active body parts & joints. Easily removable, due to its special layer on the absorbent pad, transparent wound dressing does not stick on the wound area.



- ◆ Waterproof
- ◆ STERILE unless opened or damaged & Hypoallergenic
- ◆ Transparent Elastic Polyurethane carrier
- ◆ Long-term powerful adhesion & Air permeable
- ◆ Hygienic. All sterile dressings are packed individually
- ◆ Non-absorbent wound pad
- ◆ Suitable for x-ray
- ◆ Single use

Code	Size	Quantity / Box	Quantity / Carton
21507	5cm x 7,5cm	50 pcs	1200 pcs
21509	5cm x 9cm	50 pcs	800 pcs
21910	9cm x 10cm	25 pcs	500 pcs
21915	9cm x 15cm	25 pcs	400 pcs
21920	9cm x 20cm	25 pcs	350 pcs
21925	9cm x 25cm	25 pcs	250 pcs
21930	9cm x 30cm	25 pcs	200 pcs

Store in a cool and dry place.