


**It has high quality transparent and ultra thin film with high elasticity.**

**Transparent structure permits skin to breathe and prevents gathering of excessive moisture under the film and decreases infection possibility and provides an acceptable environment around the cuttings.**

**Blood and similar liquids are suitable environments for bacteria.**

**HI - TEC Surgical drape covers the operation site with its adhesive area and does not permit to effect operation area by bacteria, forms a barrier against bacteria and does not permit the transfer of bacteria. If the operation takes longer than 4-6 hours, it is recommended to use a secondary drape as the adhesive efficiency decreases by time. Surgical Incise Drape is subjected to ethylene oxide sterilization.**


CAT. NO	ADHESIVE SIZE	FULL SIZE
D1528	15X28 cm	21X28 cm
D2020	20X20 cm	26X20 cm
D3040	30X40 cm	36X40 cm
D3535	35X35 cm	41X35 cm
D4545	45X45 cm	51X45 cm
D5645	56X45 cm	62X45 cm
D6045	60X45 cm	66X45 cm
D6060	60X60 cm	66X60 cm
D6090	60X90 cm	66X90 cm
CUSTOMIZED SIZES ARE AVAILABLE		



HI-TEC Surgical Incision Drape is used in surgical operations to prevent the spreading of skin bacteria and post-op operations.

After several clinical studies show that the skin bacterias have big role of post-op infections.

HI-TEC surgical drape is used in all type of surgeries which caused by post-op skin bacterias.



It is made of cuttable polyurethane and/or polyethylene film.

Surgical Incise Drape is made of release paper and 1 layer of incise film.

Impermeable to bacteria and water

Hypoallergenic

Proper viscosity and high elasticity

Bone, Abdomen operation, Chest and

Heart surgical operation, Plastic operation,

Ophthalmology and Pediatrics operations are some application areas.



Incise Drape